




Data sheet information

1. PRODUCT IDENTIFICATION

Product code	Product name	Photo	Package
P2265	Blue detectable plasters mix (3 sizes)		100 pcs

2. PROPERTIES

Components	Material	Dimensions	Woundpad
Blue detectable plasters – 40 psc	PE (polietylen) Colour: Blue Scent: Odorless + aluminium film under wound pad		The aluminium woundpad is made of 0.025mm thick aluminum sheet (approximately 25x11mm) which is attached with an hypoallergenic glue to a non-adhesive woundpad. The hypoallergenic glue doesn't contain latex. And the film doesn't contain SVHC
Blue detectable plasters – 40 psc			The aluminium woundpad is made of 0.025mm thick aluminum sheet (approximately 25x15mm) which is attached with an hypoallergenic glue to a non-adhesive woundpad. The hypoallergenic glue doesn't contain latex. And the film doesn't contain SVHC
Blue detectable plasters – 20 psc			The aluminium woundpad is made of 0.025mm thick aluminum sheet (approximately 35x25mm) which is attached with an hypoallergenic glue to a non-adhesive woundpad. The hypoallergenic glue doesn't contain latex. And the film doesn't contain SVHC

3. OTHER INFORMATION

a) Used

Store in sealed rooms; Keep out of reach of any fire, glow, heat or ignition source; Maximum storage temperature max + 25°C, min -5° C; under moderate air humidity conditions, to protect the package against sunlight, water and humidity the intrusion.
Keep out of reach of oxidating products.
Meet the necessary regulations for sterile storage requirements.

b) Risk

Not harmful under normal industrial circumstances & personal hygien.
Skin irritation due to dehydration is unlikely because of the perforated & non adhesive woundpad.
Nevertheless skin irritation cannot be ruled out when used repeatedly or under long exposure.

4. REGULATORY INFORMATION

Classification: Class I Medical devices

Traceability module: Lot & expiry date on backside of the box. With possible barcode.

Regulation (EC) nr. 1935/2004 This product complies with the regulation in order to be exposed to human consumption goods.

5. DETECTABILITY

PROHACCP Sp. z o.o. informs that our metal detectable products have magnetic characteristics and can be detected by standard metal detection systems. However due to settings of different metal detection systems, PROHACCP Sp. z o.o. cannot guarantee that our products will be detected by your metal detectors. Therefore we cannot take responsibility for any non-detection of our products or their parts. Detectability of our products depends on metal detector settings, and type of food products produced (product composition, size, humidity etc.). In order to make sure that our metal detectable products will be detected we strongly recommend that you test it on your metal detectors before you use it.

Kielce, November 10, 2020

Certified on behalf of PROHACCP Sp. z o. o.

Edyta Kowalczyk

Quality Manager